

Peak Forearm Crutches, Plastic Cuff

Features:

- Light Weight Aluminium Frame
- Plastic Cuff
- Handgrip Attached to Frame by Rivet
- 3 Inches for Elbow Height Adjustment
- Locking Anti-Rattle Device on Both
Top and Bottom Sections
- Double Safety Leg Adjustment
Buttons



Product Code	Max User Height	Min User Height	Weight	Capacity
9346376062114	5' 4"	4' 7"	0.62kg	130kg
9346376062121	6' 1"	5' 4"	0.62kg	130kg
9346376062138	6' 3"	5' 6"	0.62kg	130kg