

# Brado Medium Rollator - Ivory

Product Code: C4500-MIA-FX













## Features:

- ▶ One Hand Folding Mechanism
- ▶ Secure Brake Design
- ▶ Mesh Bag Included
- ▶ Cane Holder Included
- ▶ Foldable

## ARTG Listed:

- ▶ 227600



Length (cm)		66	Folded Width (cm)		22	Front Castor Diameter (cm)		19.4/8"
Width (cm)		60	Seat Depth (cm)		19	Rear Castor Diameter (cm)		19.4/8"
Max Height (cm)		102	Seat Width (cm)		41	Max Load (kg)		150
Minimum Height (cm)		77	Seat Height from the Ground (cm)		62.5	Net Weight (kg)		7.2