













Knee Cruiser

Product Code: CTKC-DC-FX

Features:

- Ideal for Foot or Ankle Injuries
- Stable & Maneuverable
- Locking Handbrakes
- Padded & Contoured Platform
- Basket Included
- Anatomical Handles



Length (cm)		N/A	Folded Width (cm)		45	Front Castor Diameter (cm)		20.3/8"
Width (cm)		85	Seat Depth (cm)		N/A	Rear Castor Diameter (cm)		20.3/8"
Max Height (cm)		78	Max Seat Height from the Ground (cm)		57	Max Load (kg)		135
Minimum Height (cm)		99	Min Seat Height from the Ground (cm)		43	Net Weight (kg)		11.5